

产品名称: HoxC6 Rabbit Monoclonal Antibody
产品货号: AMRe84117



产品概述 (Summary)

产品名称 (Production Name)	HoxC6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HoxC6
别名 (Alternative Names)	CP25 ; HOX3; HOX3C; Hoxc6;; HOX C6
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P09630.A synthesized peptide derived from human HOX C6

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 27 kDa ; Observed MW: 26 kDa

研究背景 (Background)

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific


产品名称: HoxC6 Rabbit Monoclonal Antibody
产品货号: AMRe84117




positional identities on the anterior-posterior axis.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




Western blot analysis of HoxC6 expression in MCF7 cell lysate.


产品名称: HoxC6 Rabbit Monoclonal Antibody
产品货号: AMRe84117




Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:50 dilution.




Immunofluorescent analysis using the Antibody at 1:150 dilution.



产品名称: HoxC6 Rabbit Monoclonal Antibody
产品货号: AMRe84117



Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:150 dilution.

注意事项 (Note)

For research use only .