

产品名称: HSC70 Interacting Protein HIP Rabbit Monoclonal Antibody
产品货号: AMRe84123



产品概述 (Summary)

产品名称 (Production Name)	HSC70 Interacting Protein HIP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSC70 Interacting Protein HIP
别名 (Alternative Names)	HIP; HOP; P48; AAG2; SNC6; HSPABP; FAM10A1; FAM10A4; HSPABP1; ST13;;Hip
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P50502.A synthesized peptide derived from human Hip

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 41 kDa ; Observed MW: 47 kDa

产品名称: HSC70 Interacting Protein HIP Rabbit Monoclonal Antibody
产品货号: AMRe84123




研究背景 (Background)


One HIP oligomer binds the ATPase domains of at least two HSC70 molecules dependent on activation of the HSC70 ATPase by HSP40. Stabilizes the ADP state of HSC70 that has a high affinity for substrate protein.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of HSC70 Interacting Protein HIP expression in 293T cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: HSC70 Interacting Protein HIP Rabbit Monoclonal
Antibody
产品货号: AMRe84123



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .