

产品名称: MAG Rabbit Monoclonal Antibody
产品货号: AMRe84125



产品概述 (Summary)

产品名称 (Production Name)	MAG Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAG
别名 (Alternative Names)	GMA; MAG; S MAG; SIGLEC4A; SPG75;;MAG
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P20916.A synthesized peptide derived from human MAG

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 69 kDa ; Observed MW: 100 kDa

研究背景 (Background)

Adhesion molecule in postnatal neural development that mediates sialic-acid dependent cell-cell interactions between


产品名称: MAG Rabbit Monoclonal Antibody
产品货号: AMRe84125




neuronal and myelinating cells. Preferentially binds to alpha-2,3-linked sialic acid.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MAG expression in Rat brain lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .