

产品名称: Gemin2 Rabbit Monoclonal Antibody
产品货号: AMRe84127



产品概述 (Summary)

产品名称 (Production Name)	Gemin2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Gemin2
别名 (Alternative Names)	gemin2; SIP 1; SIP1; SIP1 delta; SMN interacting protein 1; Survival interacting protein 1;;Gemin 2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O14893.A synthesized peptide derived from human Gemin 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 32 kDa ; Observed MW: 31 kDa

研究背景 (Background)


产品名称: Gemin2 Rabbit Monoclonal Antibody
产品货号: AMRe84127




The SMN complex plays an essential role in spliceosomal snRNP assembly in the cytoplasm and is required for pre-mRNA splicing in the nucleus.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Gemin 2 expression in HepG2 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:3K dilution)

产品名称: Gemin2 Rabbit Monoclonal Antibody
产品货号: AMRe84127



注意事项 (Note)

For research use only .