

产品名称: NR1D1 Rabbit Monoclonal Antibody
产品货号: AMRe84133



产品概述 (Summary)

产品名称 (Production Name)	NR1D1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NR1D1
别名 (Alternative Names)	EAR1; hRev; Nr1d1; Rev erbAalpha; Reverba; THRA1; THRAL;;NR1D1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P20393.A synthesized peptide derived from human NR1D1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 67 kDa ; Observed MW: 58 kDa

研究背景 (Background)


产品名称: NR1D1 Rabbit Monoclonal Antibody
产品货号: AMRe84133



Transcriptional repressor which coordinates circadian rhythm and metabolic pathways in a heme-dependent manner. Integral component of the complex transcription machinery that governs circadian rhythmicity and forms a critical negative limb of the circadian clock by directly repressing the expression of core clock components BMAL1, CLOCK and CRY1.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NR1D1 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .