

产品名称: RPS19 Rabbit Monoclonal Antibody
产品货号: AMRe84134



产品概述 (Summary)

产品名称 (Production Name)	RPS19 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPS19
别名 (Alternative Names)	DBA; DBA1; Ribosomal protein S19; rps19; S19;;RPS19
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P39019.A synthesized peptide derived from human RPS19

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	16 kDa

研究背景 (Background)

Component of the small ribosomal subunit.The ribosome is a large ribonucleoprotein complex responsible for the synthesis


产品名称: RPS19 Rabbit Monoclonal Antibody
产品货号: AMRe84134



of proteins in the cell. Required for pre-rRNA processing and maturation of 40S ribosomal subunits.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RPS19 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .