

产品名称: MAT2A Rabbit Monoclonal Antibody
产品货号: AMRe84135



产品概述 (Summary)

产品名称 (Production Name)	MAT2A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAT2A
别名 (Alternative Names)	AMS2; MAT 2; MAT 2A; MAT II; Mat2a; MATA2; MATII; SAMS 2; SAMS2;;MAT2A
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P31153.A synthesized peptide derived from human MAT2A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 44 kDa ; Observed MW: 48 kDa

研究背景 (Background)


Catalyzes the formation of S-adenosylmethionine from methionine and ATP.

产品名称: MAT2A Rabbit Monoclonal Antibody
产品货号: AMRe84135




研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of MAT2A expression in (1) Ramos cell lysate;(2) Mouse spleen lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

产品名称: MAT2A Rabbit Monoclonal Antibody
产品货号: AMRe84135



For research use only .