

产品名称: GABA A Receptor beta 2 Rabbit Monoclonal Antibody
产品货号: AMRe84136



产品概述 (Summary)

产品名称 (Production Name)	GABA A Receptor beta 2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GABA A Receptor beta 2
别名 (Alternative Names)	GABA; GABA(A) receptor beta 2; GABA(A) receptor subunit beta-2; GABRB2;; GABRB2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P47870.A synthesized peptide derived from human GABRB2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	59 kDa

产品名称: GABA A Receptor beta 2 Rabbit Monoclonal Antibody
产品货号: AMRe84136




研究背景 (Background)

Beta subunit of the heteropentameric ligand-gated chloride channel gated by gamma-aminobutyric acid (GABA), a major inhibitory neurotransmitter in the brain. GABA-gated chloride channels, also named GABA(A) receptors (GABAAR), consist of five subunits arranged around a central pore and contain GABA active binding site(s) located at the alpha and beta subunit interface(s).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GABA A Receptor beta 2 expression in A549 cell lysate.

注意事项 (Note)

For research use only .