

产品名称: ILF2 Rabbit Monoclonal Antibody
产品货号: AMRe84140



产品概述 (Summary)

产品名称 (Production Name)	ILF2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ILF2
别名 (Alternative Names)	ILF 2; MGC8391; NF45;;ILF2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q12905.A synthesized peptide derived from human ILF2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	43 kDa

研究背景 (Background)

Appears to function predominantly as a heterodimeric complex with ILF3. This complex may regulate transcription of the


产品名称: ILF2 Rabbit Monoclonal Antibody
产品货号: AMRe84140



IL2 gene during T-cell activation. It can also promote the formation of stable DNA-dependent protein kinase holoenzyme complexes on DNA.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ILF2 expression in A431 cell lysate.

注意事项 (Note)

For research use only .