

产品名称: MAGOH Rabbit Monoclonal Antibody
产品货号: AMRe84142



产品概述 (Summary)

产品名称 (Production Name)	MAGOH Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAGOH
别名 (Alternative Names)	magoh; MAGOHA; Protein mago nashi homolog;;MAGOH
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P61326.A synthesized peptide derived from human MAGOH

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	17 kDa

研究背景 (Background)

Component of a splicing-dependent multiprotein exon junction complex (EJC) deposited at splice junction on mRNAs. The


产品名称: MAGOH Rabbit Monoclonal Antibody
产品货号: AMRe84142



EJC is a dynamic structure consisting of a few core proteins and several more peripheral nuclear and cytoplasmic associated factors that join the complex only transiently either during EJC assembly or during subsequent mRNA metabolism.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MAGOH expression in HeLa cell lysate.

注意事项 (Note)

For research use only .