

产品名称: Nogo Receptor Rabbit Monoclonal Antibody
产品货号: AMRe84151



产品概述 (Summary)

产品名称 (Production Name)	Nogo Receptor Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Nogo Receptor
别名 (Alternative Names)	NgR; Nogo 66 receptor; Nogo receptor; Nogor; Reticulon 4 receptor; Rtn4r;;Reticulon 4 receptor
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9BZR6.A synthesized peptide derived from human Reticulon 4 receptor

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 51 kDa ; Observed MW: 66 kDa

研究背景 (Background)


产品名称: Nogo Receptor Rabbit Monoclonal Antibody
产品货号: AMRe84151



Receptor for RTN4, OMG and MAG, Functions as a receptor for the sialylated gangliosides GT1b and GM1. Besides, functions as a receptor for chondroitin sulfate proteoglycans.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Nogo Receptor expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .