

产品名称: E2F6 Rabbit Monoclonal Antibody
产品货号: AMRe84152



产品概述 (Summary)

产品名称 (Production Name)	E2F6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	E2F6
别名 (Alternative Names)	E2F 6; E2F-6; E2F6; EMA;;E2F6
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O75461.A synthesized peptide derived from human E2F6

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	32 kDa

研究背景 (Background)

Inhibitor of E2F-dependent transcription. Binds DNA cooperatively with DP proteins through the E2 recognition site, 5'-


产品名称: E2F6 Rabbit Monoclonal Antibody
产品货号: AMRe84152



TTTC[CG]CGC-3'. Has a preference for the 5'-TTTCCCGC-3' E2F recognition site. E2F-6 lacks the transcriptional activation and pocket protein binding domains.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of E2F6 expression in K562 cell lysate.

注意事项 (Note)

For research use only .