

产品名称: KLRG1 Rabbit Monoclonal Antibody
产品货号: AMRe84154



产品概述 (Summary)

产品名称 (Production Name)	KLRG1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KLRG1
别名 (Alternative Names)	2F1 Ag; CLEC15A; KLRG1; MAFA 2F1; MAFA L; MAFA like; MAFAL;;KLRG1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q96E93.A synthesized peptide derived from human KLRG1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 22 kDa ; Observed MW: 38 kDa

研究背景 (Background)

Plays an inhibitory role on natural killer (NK) cells and T-cell functions upon binding to their non-MHC ligands. May


产品名称: KLRG1 Rabbit Monoclonal Antibody
产品货号: AMRe84154




mediate missing self recognition by binding to a highly conserved site on classical cadherins, enabling it to monitor expression of E-cadherin/CDH1, N-cadherin/CDH2 and R-cadherin/CDH4 on target cells.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of KLRG1 expression in Jurkat cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: KLRG1 Rabbit Monoclonal Antibody
产品货号: AMRe84154



注意事项 (Note)

For research use only .