

产品名称: AKR1C1/AKR1C2 Rabbit Monoclonal Antibody
产品货号: AMRe84161



产品概述 (Summary)

产品名称 (Production Name)	AKR1C1/AKR1C2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AKR1C1/AKR1C2
别名 (Alternative Names)	AKR1C1; AKR1C2;;AKR1C1/2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q04828/P52895.A synthesized peptide derived from human AKR1C1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	37 kDa

研究背景 (Background)

In the liver and intestine, may have a role in the transport of bile. May have a role in monitoring the intrahepatic bile acid


产品名称: AKR1C1/AKR1C2 Rabbit Monoclonal Antibody
产品货号: AMRe84161




concentration. Has a low bile-binding ability. May play a role in myelin formation.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of AKR1C1/AKR1C2 expression in HepG2 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: AKR1C1/AKR1C2 Rabbit Monoclonal Antibody
产品货号: AMRe84161



注意事项 (Note)

For research use only .