

产品名称: PDE4B Rabbit Monoclonal Antibody
产品货号: AMRe84168



产品概述 (Summary)

产品名称 (Production Name)	PDE4B Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PDE4B
别名 (Alternative Names)	DPDE4; PDE 4B; PDE IVB; PDE32; Pde4b; PDE4B5; PDEIVB; Phosphodiesterase 4B;;PDE4B
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q07343.A synthesized peptide derived from human PDE4B

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 83 kDa ; Observed MW: 95 kDa

研究背景 (Background)


产品名称: PDE4B Rabbit Monoclonal Antibody
产品货号: AMRe84168



Hydrolyzes the second messenger cAMP, which is a key regulator of many important physiological processes. May be involved in mediating central nervous system effects of therapeutic agents ranging from antidepressants to antiasthmatic and anti-inflammatory agents.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PDE4B expression in A549 cell lysate.

注意事项 (Note)

For research use only .