

产品名称: CYP26A1 Rabbit Monoclonal Antibody
产品货号: AMRe84170



产品概述 (Summary)

产品名称 (Production Name)	CYP26A1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CYP26A1
别名 (Alternative Names)	CP26; CYP26; cyp26a1; hP450RAI; P450RAI; P450RAI1;;Cytochrome P450 26A1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O43174.A synthesized peptide derived from human Cytochrome P450 26A1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 56 kDa ; Observed MW: 50 kDa

研究背景 (Background)

A cytochrome P450 monooxygenase involved in the metabolism of retinoates (RAs), the active metabolites of vitamin A,


产品名称: CYP26A1 Rabbit Monoclonal Antibody
产品货号: AMRe84170




and critical signaling molecules in animals.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of CYP26A1 expression in HepG2 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: CYP26A1 Rabbit Monoclonal Antibody
产品货号: AMRe84170



注意事项 (Note)

For research use only .