

产品名称: ATP1A2 Rabbit Monoclonal Antibody
产品货号: AMRe84171



产品概述 (Summary)

产品名称 (Production Name)	ATP1A2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATP1A2
别名 (Alternative Names)	Atp1a2; FHM2; KIAA0778; MHP2;;ATP1A2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P50993.A synthesized peptide derived from human ATP1A2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 112 kDa ; Observed MW: 102 kDa

研究背景 (Background)

This is the catalytic component of the active enzyme, which catalyzes the hydrolysis of ATP coupled with the exchange of


产品名称: ATP1A2 Rabbit Monoclonal Antibody
产品货号: AMRe84171



sodium and potassium ions across the plasma membrane. This action creates the electrochemical gradient of sodium and potassium, providing the energy for active transport of various nutrients.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ATP1A2 expression in Human fetal heart lysate.

注意事项 (Note)

For research use only .