

产品名称: Septin2 Rabbit Monoclonal Antibody
产品货号: AMRe84176



产品概述 (Summary)

产品名称 (Production Name)	Septin2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Septin2
别名 (Alternative Names)	DIFF6; hNEDD5; KIAA0158; NEDD5; Pnutl3; SEPT2; septin 2;;Septin 2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q15019.A synthesized peptide derived from human Septin 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	41 kDa

研究背景 (Background)


产品名称: Septin2 Rabbit Monoclonal Antibody
产品货号: AMRe84176



Filament-forming cytoskeletal GTPase. Forms a filamentous structure with SEPTIN12, SEPTIN6, SEPTIN2 and probably SEPTIN4 at the sperm annulus which is required for the structural integrity and motility of the sperm tail during postmeiotic differentiation

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Septin 2 expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .