

产品名称: Sprouty4 Rabbit Monoclonal Antibody
产品货号: AMRe84177



产品概述 (Summary)

产品名称 (Production Name)	Sprouty4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Sprouty4
别名 (Alternative Names)	Sprouty homolog 4; SPRY 4;; Sprouty 4
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9C004.A synthesized peptide derived from human Sprouty 4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	33 kDa

研究背景 (Background)

Suppresses the insulin receptor and EGFR-transduced MAPK signaling pathway, but does not inhibit MAPK activation by a


产品名称: Sprouty4 Rabbit Monoclonal Antibody
产品货号: AMRe84177



constitutively active mutant Ras. Probably impairs the formation of GTP-Ras. Inhibits Ras-independent, but not Ras-dependent, activation of RAF1.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Sprouty 4 expression in MCF7 cell lysate.

注意事项 (Note)

For research use only .