

产品名称: DDX6 Rabbit Monoclonal Antibody
产品货号: AMRe84178



产品概述 (Summary)

产品名称 (Production Name)	DDX6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DDX6
别名 (Alternative Names)	DDX6; HLR2; P54; RCK;;DDX6
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P26196.A synthesized peptide derived from human DDX6

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	54 kDa

研究背景 (Background)

DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are


产品名称: DDX6 Rabbit Monoclonal Antibody
产品货号: AMRe84178



implicated in a number of cellular processes involving alteration of RNA secondary structure such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of DDX6 expression in K562 cell lysate.

注意事项 (Note)

For research use only .