

产品名称: GABA A Receptor alpha 5 Rabbit Monoclonal Antibody
产品货号: AMRe84183



产品概述 (Summary)

产品名称 (Production Name)	GABA A Receptor alpha 5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GABA A Receptor alpha 5
别名 (Alternative Names)	GAA 5; GAA5; GABRA 5; Gabra5;; GABRA5
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P31644.A synthesized peptide derived from human GABRA5

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 52 kDa ; Observed MW: 70 kDa

研究背景 (Background)


产品名称: GABA A Receptor alpha 5 Rabbit Monoclonal Antibody
产品货号: AMRe84183




GABA, the major inhibitory neurotransmitter in the vertebrate brain, mediates neuronal inhibition by binding to the GABA/benzodiazepine receptor and opening an integral chloride channel.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of GABA A Receptor alpha 5 expression in 293T cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: GABA A Receptor alpha 5 Rabbit Monoclonal
Antibody
产品货号: AMRe84183



注意事项 (Note)

For research use only .