

产品名称: Choline Acetyltransferase Rabbit Monoclonal Antibody
产品货号: AMRe84185



产品概述 (Summary)

产品名称 (Production Name)	Choline Acetyltransferase Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Choline Acetyltransferase
别名 (Alternative Names)	Acetyl CoA choline O; ChAT; CHOACTase; Choline acetylase; choline acetyltransferase; Choline O acetyltransferase; Choline O-acetyltransferase; CMS1A; CMS1A2;;Choline O cetyltransferase
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q03059(mouse)P32738(Rat)(guineapig:P13222)(human:P28329).A synthesized peptide derived from mouse Choline O cetyltransferase

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
----------------------	----------------------------------

产品名称: Choline Acetyltransferase Rabbit Monoclonal Antibody
产品货号: AMRe84185




蛋白分子量 (Molecular Weight) 72 kDa

研究背景 (Background)


Catalyzes the reversible synthesis of acetylcholine (ACh) from acetyl CoA and choline at cholinergic synapses.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of Choline Acetyltransferase expression in Mouse brain cell lysate.

注意事项 (Note)

For research use only .