

产品名称: PKC delta+PKCtheta Rabbit Monoclonal Antibody
产品货号: AMRe84187



产品概述 (Summary)

产品名称 (Production Name)	PKC delta+PKCtheta Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PKC delta+PKCtheta ALPS3; CVID9; MAY1; nPKC-delta; nPKC-theta; PKCD; PRKCD; PRKCQ; PRKCT;
别名 (Alternative Names)	Protein kinase C delta type; Protein kinase C theta type; SDK1; Sphingosine-dependent protein kinase-1; Tyrosine-protein kinase PRKCD;;PKC delta/theta
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q05655/Q04759/P28867/.A synthesized peptide derived from human PKC delta

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 78,82 kDa ; Observed MW: 78 kDa

产品名称: PKC delta+PKCtheta Rabbit Monoclonal Antibody
产品货号: AMRe84187




研究背景 (Background)


Calcium-independent, phospholipid- and diacylglycerol (DAG)-dependent serine/threonine-protein kinase that plays contrasting roles in cell death and cell survival by functioning as a pro-apoptotic protein during DNA damage-induced apoptosis, but acting as an anti-apoptotic protein during cytokine receptor-initiated cell death, is involved in tumor suppression as well as survival of several cancers, is required for oxygen radical production by NADPH oxidase and acts as positive or negative regulator in platelet functional responses.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of PKC delta + PKC theta expression in (1) HeLa cell lysate; (2) Raw264.7 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: PKC delta+PKCtheta Rabbit Monoclonal Antibody
产品货号: AMRe84187



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .