

产品名称: PAR6 Rabbit Monoclonal Antibody
产品货号: AMRe84188



产品概述 (Summary)

产品名称 (Production Name)	PAR6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PAR6
别名 (Alternative Names)	PAR6; PAR6C; TAX40; PAR-6A; TIP-40; PAR6alpha;;PAR6 alpha
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9NPB6.A synthesized peptide derived from human PAR6 alpha

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	37 kDa

研究背景 (Background)

Adapter protein involved in asymmetrical cell division and cell polarization processes. Probably involved in the formation of


产品名称: PAR6 Rabbit Monoclonal Antibody
产品货号: AMRe84188



epithelial tight junctions. Association with PARD3 may prevent the interaction of PARD3 with F11R/JAM1, thereby preventing tight junction assembly. The PARD6-PARD3 complex links GTP-bound Rho small GTPases to atypical protein kinase C proteins.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PAR6 expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .