

产品名称: SMARCC1/BAF155 Rabbit Monoclonal Antibody
产品货号: AMRe84189



产品概述 (Summary)

产品名称 (Production Name)	SMARCC1/BAF155 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC, IP
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SMARCC1/BAF155
别名 (Alternative Names)	BAF155; CRACC1; Rsc8; SMARCC1; SRG3; SWI3;; SMARCC1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q92922.A synthesized peptide derived from human SMARCC1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 123 kDa ; Observed MW: 155 kDa

研究背景 (Background)

Involved in transcriptional activation and repression of select genes by chromatin remodeling (alteration of DNA-


产品名称: SMARCC1/BAF155 Rabbit Monoclonal Antibody
产品货号: AMRe84189



nucleosome topology). Component of SWI/SNF chromatin remodeling complexes that carry out key enzymatic activities, changing chromatin structure by altering DNA-histone contacts within a nucleosome in an ATP-dependent manner. May stimulate the ATPase activity of the catalytic subunit of the complex.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of SMARCC1 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .