

产品名称: GRASP65 Rabbit Monoclonal Antibody
产品货号: AMRe84190



产品概述 (Summary)

产品名称 (Production Name)	GRASP65 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRASP65
别名 (Alternative Names)	GOLPH5; GORASP1; GRASP65;;GORASP1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9BQQ3.A synthesized peptide derived from human GORASP1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 46 kDa ; Observed MW: 65 kDa

研究背景 (Background)


产品名称: GRASP65 Rabbit Monoclonal Antibody
产品货号: AMRe84190




Key structural protein of the Golgi apparatus. The membrane cisternae of the Golgi apparatus adhere to each other to form stacks, which are aligned side by side to form the Golgi ribbon. Acting in concert with GORASP2/GRASP55, is required for the formation and maintenance of the Golgi ribbon, and may be dispensable for the formation of stacks.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GRASP65 expression in MCF7 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .