

产品名称: TXNRD2 Rabbit Monoclonal Antibody
产品货号: AMRe84192



产品概述 (Summary)

产品名称 (Production Name)	TXNRD2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TXNRD2
别名 (Alternative Names)	SelZ; TR 3; TR; TR beta; TR3; TRXR2; TXNRD 2; Txnrd2;;TXNRD2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9NNW7.A synthesized peptide derived from human TXNRD2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	57 kDa

研究背景 (Background)

Involved in the control of reactive oxygen species levels and the regulation of mitochondrial redox homeostasis.Maintains


产品名称: TXNRD2 Rabbit Monoclonal Antibody
产品货号: AMRe84192




thioredoxin in a reduced state. May play a role in redox-regulated cell signaling.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TXNRD2 expression in LnCap cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .