

产品名称: ERp29 Rabbit Monoclonal Antibody
产品货号: AMRe84193



产品概述 (Summary)

产品名称 (Production Name)	ERp29 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ERp29
别名 (Alternative Names)	ERp28; ERP29_HUMAN; ERp31; HEL S 107; PDI-DB; PDIA9;;ERp29
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P30040.A synthesized peptide derived from human ERp29

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	29 kDa

研究背景 (Background)

Proper protein folding and post-translational modifications are essential for secretory protein export out of the


产品名称: ERp29 Rabbit Monoclonal Antibody
产品货号: AMRe84193



endoplasmic reticulum. This task is accomplished by chaperone proteins such as protein disulfide isomerase (PDI), GRP94, and BiP.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ERp29 expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .