

产品名称: ATAD2 Rabbit Monoclonal Antibody
产品货号: AMRe84196



产品概述 (Summary)

产品名称 (Production Name)	ATAD2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATAD2
别名 (Alternative Names)	ANCCA; Atad2; CT137; L16; PRO2000;;ATAD2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q6PL18.A synthesized peptide derived from human ATAD2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 159 kDa ; Observed MW: 180 kDa

研究背景 (Background)

May be a transcriptional coactivator of the nuclear receptor ESR1 required to induce the expression of a subset of estradiol


产品名称: ATAD2 Rabbit Monoclonal Antibody
产品货号: AMRe84196




target genes, such as CCND1, MYC and E2F1. May play a role in the recruitment or occupancy of CREBBP at some ESR1 target gene promoters. May be required for histone hyperacetylation. Involved in the estrogen-induced cell proliferation and cell cycle progression of breast cancer cells.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




Western blot analysis of ATAD2 expression in MCF-7 cell lysate.



Immunofluorescent analysis using the Antibody at 1:150 dilution.

产品名称: ATAD2 Rabbit Monoclonal Antibody
产品货号: AMRe84196



Immunofluorescent analysis using the Antibody at 1:50 dilution.

注意事项 (Note)

For research use only .