

产品名称: PARP2 Rabbit Monoclonal Antibody
产品货号: AMRe84198



产品概述 (Summary)

产品名称 (Production Name)	PARP2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PARP2
别名 (Alternative Names)	pADPRT-2; hPARP-2; ADPRT2; ADPRTL2;; PARP 2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UGN5.A synthesized peptide derived from human PARP 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	66 kDa

研究背景 (Background)

Poly-ADP-ribosyltransferase that mediates poly-ADP-ribosylation of proteins and plays a key role in DNA repair. Mediates


产品名称: PARP2 Rabbit Monoclonal Antibody
产品货号: AMRe84198




glutamate, aspartate or serine ADP-ribosylation of proteins: the ADP-D-ribosyl group of NAD⁺ is transferred to the acceptor carboxyl group of target residues and further ADP-ribosyl groups are transferred to the 2'-position of the terminal adenosine moiety, building up a polymer with an average chain length of 20-30 units.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PARP2 expression in Raji cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .