

产品名称: ELMO1 Rabbit Monoclonal Antibody
产品货号: AMRe84199



产品概述 (Summary)

产品名称 (Production Name)	ELMO1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ELMO1
别名 (Alternative Names)	CED12; Ced12 homolog 1; ELMO 1;; ELMO1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q92556.A synthesized peptide derived from human ELMO1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	84 kDa

研究背景 (Background)

Involved in cytoskeletal rearrangements required for phagocytosis of apoptotic cells and cell motility. Acts in association


产品名称: ELMO1 Rabbit Monoclonal Antibody
产品货号: AMRe84199



with DOCK1 and CRK. Was initially proposed to be required in complex with DOCK1 to activate Rac Rho small GTPases. May enhance the guanine nucleotide exchange factor (GEF) activity of DOCK1.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ELMO1 expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .