

产品名称: MYL9 Rabbit Monoclonal Antibody
产品货号: AMRe84200



产品概述 (Summary)

产品名称 (Production Name)	MYL9 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MYL9
别名 (Alternative Names)	MLC2; CMH10;;Myosin light chain 9
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P24844.A synthesized peptide derived from human Myosin light chain 9

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 20 kDa ; Observed MW: 18 kDa

研究背景 (Background)

Myosin regulatory subunit that plays an important role in regulation of both smooth muscle and nonmuscle cell contractile


产品名称: MYL9 Rabbit Monoclonal Antibody
产品货号: AMRe84200




activity via its phosphorylation. Implicated in cytokinesis, receptor capping, and cell locomotion.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:1k dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1k dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1k dilution for 1 hour at room temperature.

产品名称: MYL9 Rabbit Monoclonal Antibody
产品货号: AMRe84200



注意事项 (Note)

For research use only .