

产品名称: SEC23 Rabbit Monoclonal Antibody
产品货号: AMRe84205



产品概述 (Summary)

产品名称 (Production Name)	SEC23 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SEC23
别名 (Alternative Names)	CLSD; sec23a;;Sec23A
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q15436.A synthesized peptide derived from human Sec23A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	86 kDa

研究背景 (Background)

Component of the coat protein complex II (COPII) which promotes the formation of transport vesicles from the


产品名称: SEC23 Rabbit Monoclonal Antibody
产品货号: AMRe84205



endoplasmic reticulum (ER). The coat has two main functions, the physical deformation of the endoplasmic reticulum membrane into vesicles and the selection of cargo molecules for their transport to the Golgi complex. Required for the translocation of insulin-induced glucose transporter SLC2A4/GLUT4 to the cell membrane (By similarity).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of SEC23 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .