

产品名称: TRAF6BP Rabbit Monoclonal Antibody
产品货号: AMRe84207



产品概述 (Summary)

产品名称 (Production Name)	TRAF6BP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRAF6BP
别名 (Alternative Names)	CALCOCO3; D6Ertd404e; D6Ertd772e; PRO0105; T6BP; TAX1BP1; tax1bp1b; TXBP151;; Tax1 binding protein 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q86VP1.A synthesized peptide derived from human Tax1 binding protein 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	91 kDa

研究背景 (Background)


产品名称: TRAF6BP Rabbit Monoclonal Antibody
产品货号: AMRe84207




Ubiquitin-binding adapter that participates in inflammatory, antiviral and innate immune processes as well as selective autophagy regulation. Plays a key role in the negative regulation of NF-kappa-B and IRF3 signalings by acting as an adapter for the ubiquitin-editing enzyme A20/TNFAIP3 to bind and inactivate its substrates.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of TRAF6BP expression in HepG2 cell lysate.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: TRAF6BP Rabbit Monoclonal Antibody
产品货号: AMRe84207



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .