

产品名称: TCF21 Rabbit Monoclonal Antibody
产品货号: AMRe84213



产品概述 (Summary)

产品名称 (Production Name)	TCF21 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TCF21
别名 (Alternative Names)	Capsulin; Epicardin; POD 1; POD1; TCF-21; tcf21;;TCF21
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O43680.A synthesized peptide derived from human TCF21

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, FC 1:20-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 20 kDa ; Observed MW: 26 kDa

研究背景 (Background)

Involved in epithelial-mesenchymal interactions in kidney and lung morphogenesis that include epithelial differentiation


产品名称: TCF21 Rabbit Monoclonal Antibody
产品货号: AMRe84213



and branching morphogenesis. May play a role in the specification or differentiation of one or more subsets of epicardial cell types.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TCF21 expression in 293 cell lysate.

注意事项 (Note)

For research use only .