

产品名称: VPS26 Rabbit Monoclonal Antibody
产品货号: AMRe84214



产品概述 (Summary)

产品名称 (Production Name)	VPS26 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	VPS26
别名 (Alternative Names)	HB58; Hbeta58; hVPS26; PEP8A; VPS 26; VPS26A;;VPS26A
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O75436.A synthesized peptide derived from human VPS26A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	38 kDa

研究背景 (Background)

Acts as a component of the retromer cargo-selective complex (CSC). The CSC is believed to be the core functional


产品名称: VPS26 Rabbit Monoclonal Antibody
产品货号: AMRe84214



component of retromer or respective retromer complex variants acting to prevent missorting of selected transmembrane cargo proteins into the lysosomal degradation pathway.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of VPS26 expression in A431 cell lysate.

注意事项 (Note)

For research use only .