

产品名称: DDX17 Rabbit Monoclonal Antibody
产品货号: AMRe84223



产品概述 (Summary)

产品名称 (Production Name)	DDX17 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DDX17
别名 (Alternative Names)	Ddx17; P72; RH70;;DDX17
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q92841.A synthesized peptide derived from human DDX17

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 80 kDa ; Observed MW: 70,80 kDa

研究背景 (Background)


产品名称: DDX17 Rabbit Monoclonal Antibody
产品货号: AMRe84223




As an RNA helicase, unwinds RNA and alters RNA structures through ATP binding and hydrolysis. Involved in multiple cellular processes, including pre-mRNA splicing, alternative splicing, ribosomal RNA processing and miRNA processing, as well as transcription regulation. Regulates the alternative splicing of exons exhibiting specific features.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of DDX17 expression in 293T cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: DDX17 Rabbit Monoclonal Antibody
产品货号: AMRe84223



注意事项 (Note)

For research use only .