

产品名称: USP9x Rabbit Monoclonal Antibody  
产品货号: AMRe84224



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | USP9x Rabbit Monoclonal Antibody       |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IHC,ICC/IF,ICC,FC                   |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                        |

## 产品性能 (Performance)

|                     |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                |
| 修饰 (Modification)   | Unmodified                                                                                  |
| 同种型 (Isotype)       | IgG                                                                                         |
| 克隆 (Clonality)      | Monoclonal                                                                                  |
| 形式 (Form)           | Liquid                                                                                      |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.    |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                       |

## 免疫原信息 (Immunogen)

|                        |                                                                          |
|------------------------|--------------------------------------------------------------------------|
| 基因名 (Gene Name)        | USP9x                                                                    |
| 别名 (Alternative Names) | DFFRX; FAF; FafI; Fam; hFAM; MRX99; Usp9x; X chromosome; X-linked;;USP9X |
| 基因 ID (Gene ID)        |                                                                          |
| 蛋白 ID (SwissProt ID)   | Q93008.A synthesized peptide derived from human USP9X                    |

## 产品应用 (Application)

|                          |                                                                                 |
|--------------------------|---------------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100 |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 290 kDa ; Observed MW: 292 kDa                                   |

## 研究背景 (Background)


**产品名称: USP9x Rabbit Monoclonal Antibody**  
**产品货号: AMRe84224**



Deubiquitinase involved both in the processing of ubiquitin precursors and of ubiquitinated proteins. May therefore play an important regulatory role at the level of protein turnover by preventing degradation of proteins through the removal of conjugated ubiquitin.

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of USP9x expression in HeLa cell lysate.

## 注意事项 (Note)

For research use only .