

产品名称: NFIB/NF1B2 Rabbit Monoclonal Antibody
产品货号: AMRe84230



产品概述 (Summary)

产品名称 (Production Name)	NFIB/NF1B2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NFIB/NF1B2
别名 (Alternative Names)	CTF; HMGIC/NFIB; NF-I/B; NF1-B; NF1B; NF1B2; NFI-B; NFI-RED; Nfib; NFIB2; NFIB3;;NFIB type
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O00712.A synthesized peptide derived from human NFIB type

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	47 kDa

产品名称: NFIB/NF1B2 Rabbit Monoclonal Antibody
产品货号: AMRe84230




研究背景 (Background)


Transcriptional activator of GFAP, essential for proper brain development. Recognizes and binds the palindromic sequence 5'-TTGGCNNNNNGCCAA-3' present in viral and cellular promoters and in the origin of replication of adenovirus type 2. These proteins are individually capable of activating transcription and replication.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



产品名称: NFIB/NF1B2 Rabbit Monoclonal Antibody
产品货号: AMRe84230



Western blot analysis of NFIB / NF1B2 expression in NIH/3T3 cell lysate.

注意事项 (Note)

For research use only .