

产品名称: TRPC1 Rabbit Monoclonal Antibody
产品货号: AMRe84231



产品概述 (Summary)

产品名称 (Production Name)	TRPC1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRPC1
别名 (Alternative Names)	HTRP1; Mtrp1; TRP1; TRPC1; Trrp1;;TRPC1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P48995.A synthesized peptide derived from human TRPC1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 91 kDa ; Observed MW: 120 kDa

研究背景 (Background)

Thought to form a receptor-activated non-selective calcium permeant cation channel. Probably is operated by a


产品名称: TRPC1 Rabbit Monoclonal Antibody
产品货号: AMRe84231



phosphatidylinositol second messenger system activated by receptor tyrosine kinases or G-protein coupled receptors.
Seems to be also activated by intracellular calcium store depletion.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TRPC1 expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .