

产品名称: Scramblase 1 Rabbit Monoclonal Antibody
产品货号: AMRe84234



产品概述 (Summary)

产品名称 (Production Name)	Scramblase 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Scramblase 1
别名 (Alternative Names)	MmTRA1a; MmTRA1b; Nor1; PLSCR 1; Scramblase1; Tra1; Tra1a; Tra1b; Tras1; Tras2;; Phospholipid scramblase 1
基因 ID (Gene ID)	O15162.A
蛋白 ID (SwissProt ID)	synthesized peptide derived from human Phospholipid scramblase 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	35 kDa

产品名称: Scramblase 1 Rabbit Monoclonal Antibody
产品货号: AMRe84234




研究背景 (Background)

Catalyzes calcium-induced ATP-independent rapid bidirectional and non-specific movement of phospholipids (lipid scrambling or lipid flip-flop) between the inner and outer leaflet of the plasma membrane resulting in collapse of the phospholipid asymmetry which leads to phosphatidylserine externalization on the cell surface.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Scramblase 1 expression in A431 cell lysate.

注意事项 (Note)

For research use only .