

产品名称: COPA Rabbit Monoclonal Antibody
产品货号: AMRe84236



产品概述 (Summary)

产品名称 (Production Name)	COPA Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	COPA
别名 (Alternative Names)	Alpha COP; Alpha COPI; AlphaCOP; COP A; copA; HEP COP; Proxenin; Xenin;;COPA
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P53621.A synthesized peptide derived from human COPA

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	138 kDa

研究背景 (Background)


产品名称: COPA Rabbit Monoclonal Antibody
产品货号: AMRe84236



The coatomer is a cytosolic protein complex that binds to dilysine motifs and reversibly associates with Golgi non-clathrin-coated vesicles, which further mediate biosynthetic protein transport from the ER, via the Golgi up to the trans Golgi network.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of COPA expression in MCF7 cell lysate.

注意事项 (Note)

For research use only .