

产品名称: Mov10 Rabbit Monoclonal Antibody
产品货号: AMRe84239



产品概述 (Summary)

产品名称 (Production Name)	Mov10 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Mov10
别名 (Alternative Names)	gb110; MOV 10;;MOV10
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9HCE1.A synthesized peptide derived from human MOV10

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	114 kDa

研究背景 (Background)

5' to 3' RNA helicase that is involved in a number of cellular roles ranging from mRNA metabolism and translation,


产品名称: Mov10 Rabbit Monoclonal Antibody
产品货号: AMRe84239




modulation of viral infectivity, inhibition of retrotransposition, or regulation of synaptic transmission. Plays an important role in innate antiviral immunity by promoting type I interferon production. Mechanistically, specifically uses IKKepsilon/IKBKE as the mediator kinase for IRF3 activation.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: Mov10 Rabbit Monoclonal Antibody
产品货号: AMRe84239



注意事项 (Note)

For research use only .