

产品名称: PHAPI2/APRIL Rabbit Monoclonal Antibody
产品货号: AMRe84243



产品概述 (Summary)

产品名称 (Production Name)	PHAPI2/APRIL Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PHAPI2/APRIL
别名 (Alternative Names)	ANP32B; PAL31; PHAPI2; PHAPI2a; Ssp29;;ANP32B
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q92688.A synthesized peptide derived from human ANP32B

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	29 kDa

研究背景 (Background)


产品名称: PHAPI2/APRIL Rabbit Monoclonal Antibody
产品货号: AMRe84243




Multifunctional protein that is involved in the regulation of many processes including cell proliferation, apoptosis, cell cycle progression or transcription. Regulates the proliferation of neuronal stem cells, differentiation of leukemic cells and progression from G1 to S phase of the cell cycle. As negative regulator of caspase-3-dependent apoptosis, may act as an antagonist of ANP32A in regulating tissue homeostasis.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PHAPI2 / APRIL expression in (1) Jurkat cell lysate; (2) RAW264.7 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .