

产品名称: DISC1 Rabbit Monoclonal Antibody
产品货号: AMRe84244



产品概述 (Summary)

产品名称 (Production Name)	DISC1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DISC1
别名 (Alternative Names)	C1orf136; DISC1; KIAA0457; RP4-730B13.1; SCZD9;;DISC1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9NRI5.A synthesized peptide derived from human DISC1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 94 kDa ; Observed MW: 100 kDa

研究背景 (Background)


产品名称: DISC1 Rabbit Monoclonal Antibody
产品货号: AMRe84244



Involved in the regulation of multiple aspects of embryonic and adult neurogenesis. Required for neural progenitor proliferation in the ventricular/subventricular zone during embryonic brain development and in the adult dentate gyrus of the hippocampus. Participates in the Wnt-mediated neural progenitor proliferation as a positive regulator by modulating GSK3B activity and CTNNB1 abundance.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of DISC1 expression in (1) HeLa cell lysate; (2) RAW264.7 cell lysate; (3) PC-12 cell lysate.

注意事项 (Note)

For research use only .