

产品名称: NEDD4 Rabbit Monoclonal Antibody
产品货号: AMRe84246



产品概述 (Summary)

产品名称 (Production Name)	NEDD4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NEDD4
别名 (Alternative Names)	Nedd4; PIG53; RPF1;;NEDD4
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P46934.A synthesized peptide derived from human NEDD4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 149 kDa ; Observed MW: 135 kDa

研究背景 (Background)

E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and


产品名称: NEDD4 Rabbit Monoclonal Antibody
产品货号: AMRe84246




then directly transfers the ubiquitin to targeted substrates. Involved in the pathway leading to the degradation of VEGFR-2/KDR, independently of its ubiquitin-ligase activity.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NEDD4 expression in A549 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .