

产品名称: RAP1A/1B Rabbit Monoclonal Antibody
产品货号: AMRe84248



产品概述 (Summary)

产品名称 (Production Name)	RAP1A/1B Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAP1A/1B
别名 (Alternative Names)	C21KG; G 22K; K REV; KREV 1; KREV; KREV1; RAL1B; RAP 1; RAP 1A; RAP 1B; RAP1; RAP1B; SMGP21;;RAP1A/1B
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P61224/P62834.A synthesized peptide derived from human RAP1A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	21 kDa

研究背景 (Background)


产品名称: RAP1A/1B Rabbit Monoclonal Antibody
产品货号: AMRe84248



RAP1A and RAP1B belong to a family of RAS-related proteins. These proteins share approximately 50% amino acid identity with the classical RAS proteins and have numerous structural features in common.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RAP1A + RAP1B expression in HeLa cell lysate.

注意事项 (Note)

For research use only .