

产品名称: Syntrophin alpha 1 Rabbit Monoclonal Antibody
产品货号: AMRe84249



产品概述 (Summary)

产品名称 (Production Name)	Syntrophin alpha 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Syntrophin alpha 1
别名 (Alternative Names)	Alpha 1 syntrophin; LQT12; SNT1; Snta1; Syntrophin 1; TACIP1;;alpha 1 Syntrophin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q13424.A synthesized peptide derived from human alpha 1 Syntrophin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	54 kDa

产品名称: Syntrophin alpha 1 Rabbit Monoclonal Antibody
产品货号: AMRe84249




研究背景 (Background)

Adapter protein that binds to and probably organizes the subcellular localization of a variety of membrane proteins. May link various receptors to the actin cytoskeleton and the extracellular matrix via the dystrophin glycoprotein complex. Plays an important role in synapse formation and in the organization of UTRN and acetylcholine receptors at the neuromuscular synapse. Binds to phosphatidylinositol 4,5-bisphosphate.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Syntrophin alpha 1 expression in A549 cell lysate.

注意事项 (Note)

For research use only .