

产品名称: Syntaxin Rabbit Monoclonal Antibody
产品货号: AMRe84256



产品概述 (Summary)

产品名称 (Production Name)	Syntaxin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Syntaxin
别名 (Alternative Names)	EPIM; Epimorphin; EPM; HPC 1; P35 1; STX1A; Stx1b; STX1B1; STX1B2; STX2; STX2A; STX2B; STX2C; STX3; STX3A;;Syntaxin 1B
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P61266.A synthesized peptide derived from human Syntaxin 1B

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	33 kDa

研究背景 (Background)


产品名称: Syntaxin Rabbit Monoclonal Antibody
产品货号: AMRe84256




Potentially involved in docking of synaptic vesicles at presynaptic active zones. May mediate Ca^{2+} -regulation of exocytosis acrosomal reaction in sperm.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of Syntaxin expression in Human fetal brain lysate.

注意事项 (Note)

For research use only .