

产品名称: IL17A Receptor Rabbit Monoclonal Antibody
产品货号: AMRe84263



产品概述 (Summary)

产品名称 (Production Name)	IL17A Receptor Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL17A Receptor
别名 (Alternative Names)	IL-17 receptor; IL-17RA; IL17RA; interleukin-17 receptor;;IL 17 receptor A
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q96F46Human.A synthesized peptide derived from human IL 17 receptor A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 96 kDa ; Observed MW: 160 kDa

研究背景 (Background)

Receptor for IL17A and IL17F, major effector cytokines of innate and adaptive immune system involved in antimicrobial


产品名称: IL17A Receptor Rabbit Monoclonal Antibody
产品货号: AMRe84263



host defense and maintenance of tissue integrity. Receptor for IL17A. Receptor for IL17F. Binds to IL17A with higher affinity than to IL17F.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of IL17A Receptor expression in Raji cell lysate.

注意事项 (Note)

For research use only .